

Spoke Guard

Fits most wheelchair models



Spoke Guard

The Spoke Guard covers is made of impact resistant, flexible and transparent plastic with center hole 80 mm alt. 120 mm. They are suitable for most wheelchair models, even for models with brakes.



Read the user manual prior to using the product.

Care of Sweden are registered trademarks of Care of Sweden AB.
© Care of Sweden AB, 2019. All rights reserved. Reproduction prohibited.
Last updated: 2019-07-02.



MOUNTING	FACTS
<p>Clean the spokes on the wheelchair and check that the wheel is in good condition. Then place the covers on the outside of the wheel. Make sure the covers fit firmly against the spokes on the wheelchair. Attach the cover with tie strap to the pre-drilled holes. The cover shall be tightened with four attachment points. Attachment material is not included in the delivery.</p> <p>If necessary, the Spoke Guard cover can be adjusted by cutting in the material, even new holes can be drilled. Before use check that the spoke guard covers are consistent with the wheelchair's wheels.</p>	<ul style="list-style-type: none">• Material: Polyester• Packaging: 0.5% PELD 04• Cleaning: Use mild detergent and lukewarm water.