

# CuroCell UNO®

A new generation of fully automated  
alternating pressure mattress replacement systems



# CuroCell UNO®

The CuroCell UNO® is used as an aid in the prevention and treatment of pressure ulcers up to and including category 4<sup>(1)</sup>. It stands for innovation and modern design and together with quality engineering it lays the foundation for our new generation of fully automated premium mattress systems.

The pump is fully automatic with three programs – Alternating, Pulsating and Constant low pressure. The system identifies the user and manages all settings to optimize the pressure according to the user's needs. It is equipped with a large touch screen panel with a user-friendly interface.

Care of Sweden's mattresses are supplied with a four way stretch hygienic outer cover which conforms to relevant standards for hygiene<sup>(3)</sup>, quality and safety. All covers are manufactured from modern, recyclable materials and are vapour permeable<sup>(4)</sup>, which reduce the risk of skin maceration.



## OPTIONAL

- Durable bag
- Sheet holder

## SIZES

80/85/90/105/120 x 200/210/220 x 20 cm

If you need any other sizes please contact your distributor

## REFERENCES

(1) National Pressure Ulcer Advisory Panel, European Pressure Ulcer Advisory Panel and Pan Pacific Pressure Injury Alliance. Prevention and Treatment of Pressure Ulcers: Quick Reference Guide. Emily Haesler (Ed.). Cambridge Media: Osborne Park, Western Australia; 2014.

(2) "EN ISO 11201 Acoustics - Noise emitted by machinery and equipment - Determination of emission sound pressure levels", SP 2015.

(3) ISO 16603 "resistance of penetration by blood and body fluids", ISO 16604 "resistance of penetration by blood-borne pathogens", EN 14126.

(4) SS-EN ISO 15496:2004, DIN 53122-1.



Always read the instruction for use prior to usage of the product.

Care of Sweden and CuroCell UNO are trademarks of Care of Sweden AB. © Care of Sweden AB, 2016. All rights reserved. Last update: 2018-12-04.

FEATURES MATTRESS	FACTS
<ul style="list-style-type: none"> <li>• Cell on Cell construction with safety mattress</li> <li>• Heel Function</li> <li>• Static head cells</li> <li>• CPR-function for rapid deflation</li> <li>• Smart cable management</li> <li>• PVC free</li> <li>• Unique design for easy handling and usage</li> <li>• Chlorine resistant (10 %) multi stretch cover with welded seams</li> <li>• Well protected zipper</li> </ul>	<ul style="list-style-type: none"> <li>• <b>Type: Alternating pressure mattress replacement system</b></li> <li>• <b>Height: 20 cm</b></li> <li>• <b>Function: Prevention and treatment up to and including category 4<sup>(1)</sup></b></li> <li>• <b>Recommended user weight: 0-250 kg</b></li> <li>• <b>Input voltage: 100-240 V / 50-60 Hz</b></li> <li>• <b>Sound level, pump, max.: 20 dBA<sup>(2)</sup></b></li> <li>• <b>Warranty: 2 years</b></li> <li>• <b>Cleaning of cover: Wipe with cleaning agent and/or disinfectants</b></li> <li>• <b>Machine wash max 95 °C, tumble drying</b></li> </ul>
<h3>FEATURES PUMP</h3>	<div data-bbox="874 1758 986 1870" data-label="Image"> </div> <p data-bbox="1018 1778 1246 1832">Scan QR code to watch presentation video.</p> <div data-bbox="874 1899 986 2011" data-label="Image"> </div> <h3 data-bbox="1018 1890 1326 1921">MATTRESS SELECTOR</h3> <p data-bbox="1018 1924 1422 2007">- Guidelines when choosing a product See how to choose the right mattress at: <a href="http://www.careofsweden.com/mattress-selector">www.careofsweden.com/mattress-selector</a></p>
<ul style="list-style-type: none"> <li>• Fully automated pump that adapts individually</li> <li>• Large colour touch display with user-friendly interface</li> <li>• Built in instruction manual</li> <li>• 4 modes: Alternating, Pulsating, Static and Max-firm</li> <li>• Pack&amp;Go® function</li> <li>• Silent running and vibration free</li> <li>• Battery backup (optional)</li> <li>• Soft friction hangers with extra soft tip</li> <li>• Notification function for loss of air and or electricity supply</li> </ul>	